

A man with glasses, wearing a light-colored short-sleeved shirt and khaki pants, is sitting in a wooden rocking chair on a porch. He is holding an open book and reading. A large, fluffy white dog is lying on the porch floor next to him, looking towards the right. The background shows a lush green landscape with trees and a wooden barrel. The scene is bright and sunny, suggesting a peaceful outdoor setting.

LATITUDE[®]
Patient Management System

**Boston
Scientific**



LATITUDE® Patient Management System:

- The LATITUDE Communicator collects information from your device in the convenience of your home.
- The Communicator can check your cardiac device in two ways:
 - Regularly scheduled follow-ups
 - Remote monitoring between follow-ups as determined by your physician
- Your information is sent to a secure website that only your health care support team can view.
- Your physician can access your information when needed.

Your implanted device needs to be checked regularly to review information that is stored in the device and to monitor settings to ensure that the device is programmed optimally to manage your heart conditions.

The LATITUDE Patient Management system offers you and your family members greater convenience, peace of mind and the power to receive a higher standard of cardiac care.





Receiving the Technology that is right for you

There are three kinds of Communicators, touch screen wireless, push button wireless and wanded. The type of Communicator you will receive depends on the kind of cardiac device you have implanted.

The LATITUDE Patient Management system is easy to set up and use. A power cord and telephone landline is all you need to connect.

It is recommended that your LATITUDE Communicator stay near your bedside. The Communicator captures data from your device, according to a schedule set by your physician, and sends the data through your regular landline phone to a secure website where only your health care support team can view your information.



The touch screen wireless communicator

The monitoring of your heart health and specific device information will be done automatically according to a schedule set by your physician.



The push button wireless communicator

The monitoring of your heart health and specific device information will be done automatically according to a schedule set by your physician.



For patients with a wanded communicator

The wanded communicator is easy to use. The blinking action button will remind you when it's time to transmit your information according to a schedule set by your physician.

Simply hold the wand over your device, then push the blue interrogate button to transmit your information. On average, it takes 10–15 seconds to complete but may take more time if additional information is needed.

No matter which kind of Communicator you receive, you won't feel a thing when the device is being checked.

Rest assured that the Communicator only reads information from your device. It cannot change your device settings, so don't worry about pushing the wrong buttons.

Your Communicator will be shipped to your home via Federal Express in a brown Boston Scientific box



Please set up your Communicator as soon as you receive it so your clinic can immediately begin monitoring your heart health and specific device information.



Inside the box you will find:

- Your Communicator
- Phone Cord
- Power Supply
- Video Setup Guide
- Illustrated One-Page Instruction Guide
- Patient Manual

Set up is simple:

- The LATITUDE Communicator plugs into your regular phone landline and power outlet
- After plugging it in, your Communicator will guide you through the set up process

LATITUDE® Patient Support

If you or your family has questions regarding:

- Set up / Operation
- Other LATITUDE-related topics

The LATITUDE Patient Support hotline is available to you by dialing 1.866.484.3268. This live, dedicated, multi-lingual service is available to you at no additional charge.

Additional Resources for You and Your Family

Learn more about living with LATITUDE at www.aboutlatitude.com

Learn more about living with your implanted device at www.lifebeatonline.com



LATITUDE® Patient Management System from Boston Scientific CRM

The LATITUDE Patient Management system is used to remotely communicate with a compatible pulse generator device from BSC CRM and send data to a central database. The LATITUDE system is contraindicated for use with any pulse generator other than a device from BSC CRM.

The LATITUDE system is designed to tell your doctor within 24 hours if alert conditions are detected by the Communicator. However, alert notification cannot occur if:

- The Communicator is unplugged or is not able to connect to the LATITUDE system through an active phone line.
- Your device and the Communicator cannot establish and complete a communication session. This session must be initiated by you if you have a device that uses inductive telemetry (Communicator that has a wand).
- The Communicator becomes damaged or it malfunctions

Up to two weeks may go by before the LATITUDE system detects the events mentioned above, and additional time may be required for notification and resolution of the condition.

The wanded and wireless Communicator uses a radio frequency (RF) communication system to communicate with an optional weight scale and blood pressure monitor. This communication can be disrupted by electromagnetic interference. Avoid placing your Communicator next to or in the immediate vicinity of other wireless products and sources of electromagnetic energy. The wireless Communicator uses RF to also send and receive signals from the implanted device (RF enabled devices only). Using the blue Interrogate button more than as prompted by your Communicator or as instructed by your physician may lead to a decrease in the battery life of your implanted device. Your communicator is designed to be used in the continental US, Alaska, Hawaii, and Puerto Rico. These devices are available by prescription only.

(Rev. I)

**Boston
Scientific**

Delivering what's next.™

Cardiac Rhythm Management

Boston Scientific Corporation
4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA
Tel: 651.582.4000 Fax: 651.582.4166

Medical Professionals:
1.800.CARDIAC (227.3422)

Patients and Families:
1.866.484.3268

www.bostonscientific.com

Copyright © 2008 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

C7-409-1008