

Blazer Prime™ XP Temperature Ablation Catheter



Product Information


Electrode Configuration: Quadripolar

Electrode Spacing: 2.5 mm

Catheter Model No.	Shaft Size	Tip Size	Curve Style	Cable Model No.*
M004 P4500TH 0	7F	8F/8 mm (straight)	Standard	M004 651 0
M004 P4500THK2 0	7F	8F/8 mm (straight)	Large	M004 651 0
M004 P4500THN4 0	7F	8F/8 mm (straight)	Asymmetric 4	M004 651 0
M004 P4790TH 0	7F	8F/10 mm (straight)	Standard	M004 651 0
M004 P4790THK2 0	7F	8F/10 mm (straight)	Large	M004 651 0
M004 P4790THN4 0	7F	8F/10 mm (straight)	Asymmetric 4	M004 651 0
M004 P4770THK2 0	7F	8F/8 mm (VM)	Large	M004 651 0

*Note: All Blazer Prime XP Catheters are 110 cm in length; distal shaft length is 11 cm. *If an extension cable is required, use Model No. M004 613 0 (10 ft. length). All Blazer Prime XP Catheters require the use of 2 Valleylab Ground Pads (Model M004 354 0).*

Bidirectional Curve Options



Trademarks in this document are the property of Boston Scientific Corporation or its affiliates.
 © 2009 Boston Scientific Corporation or its affiliates. All rights reserved. EPT-10867_02/09 page 1 of 1

**Boston
Scientific**

For additional product information call Customer Service at (888) 272-1001
www.bostonscientific.com/electrophysiology